## **NSW HEALTH**

## Vaccine Refrigerator Protocol



## REFRIGERATOR MONITORING



**CHECK** refrigerator temperatures twice a day during operational hours and **aim for 5° C** 



**RECORD** the minimum, maximum and current temperatures on the NSW Health vaccine refrigerator temperature chart



**RESET** the thermometer after recording temperatures



**ACT** if temperatures are outside of the recommended temperature range of +2°C to +8°C or the fridge has alarmed. Follow the cold chain breach protocol steps below.



## **COLD CHAIN BREACH STEPS**

- Take corrective ACTION where possible. Is the refrigerator door closed, refrigerator plugged in/turned on. Contact engineer if a refrigerator malfunction is suspected.
- 2 ISOLATE affected vaccines/medicines and label 'DO NOT USE'. Ensure vaccines can continue to be stored between +2°C to +8°C. Vaccines may need to be transferred to an alternative monitored refrigerator or cooler.
- **REPORT breach** immediately to your manager/delegate Ph:
- **DOWNLOAD the data logger** or review the back to base or automated temperature monitoring system and investigate potential cause and duration of breach.
- NOTIFY the public health unit (1300 066 055) if temperature is <+2°C or > +8°C (excludes fluctuations up to +12°C ≤15 mins e.g. vaccine delivery). PHU will require cold chain breach reporting form, data logging and twice daily temperature readings to assess breach.
- DO NOT DISCARD vaccines until advice is received.
- OTHER MEDICATIONS Notify facility Pharmacy Service provider.

  Do NOT use medications until advice is received from Pharmacy.
- **DOCUMENT** all activity on the back of the 'NSW Vaccine Refrigerator Temperature Chart - Fortnightly,' i.e. thermostat adjustment, restocking, power outage, data logger removal etc. **SUBMIT an IIMS/ims+** for breaches resulting in vaccine wastage or patient recall and revaccination.



Order the *NSW Vaccine Refrigerator Temperature Chart - Fortnightly* from Toll Stream Direct. ORDER NUMBER: NH700227